




Interference Microscope for Razor Blade Tip Section Measurement



Microscope Service & Sales Ltd
6 The Courtyard
Furlong Road
Bourne End
Buckinghamshire SL8 5AU
Tel: +44 (0) 1628 536936
Fax: +44 (0) 1628 525296
Email: sales@microscopy.uk.com



The T10 microscope is supplied complete for T10 blade thickness measurement with:

- Mechanical stage with X movement and precise Y movement
- Blade holder
- Focusable eyepieces with graticule
- Trinocular eyepiece tube
- 50W illumination
- Interference objective and intermediate optics
- Vibration-isolation platform
- Colour camera system
- PC based image display and capture system

Microscope Service & Sales Ltd
6 The Courtyard
Furlong Road
Bourne End
Buckinghamshire SL8 5AU
Tel: +44 (0) 1628 536936
Fax: +44 (0) 1628 525296
Email: sales@microscopy.uk.com